



EMERSONTM



Branson Ultrasonic Bath

Operating Manual

Branson Ultrasonics Corp.
120 Park Ridge Road
Brookfield, CT 06804
(203) 796-0400
<http://www.bransonultrasonics.com>

BRANSON

Manual Change Information

At Branson, we strive to maintain our position as the leader in ultrasonics plastics joining, metal welding, cleaning and related technologies by continually improving our circuits and components in our equipment. These improvements are incorporated as soon as they are developed and thoroughly tested.

Information concerning any improvements will be added to the appropriate technical documentation at its next revision and printing. Therefore, when requesting service assistance for specific units, note the Revision information found on this document, and refer to the printing date which appears on this page.

Copyright and Trademark Notice

Copyright © 2023 Branson Ultrasonics Corporation. All rights reserved. Contents of this publication may not be reproduced in any form without the written permission of Branson Ultrasonics Corporation.

Other trademarks and service marks mentioned herein are held by their respective owners.

Foreword

Congratulations on your choice of a Branson Ultrasonics Corporation system!

The Branson Ultrasonic Bath system is process equipment for the cleaning of parts using ultrasonic energy. It is the newest generation of product using this sophisticated technology for a variety of customer applications. This Operating Manual is part of the documentation set for this system, and should be kept with the equipment.

Thank you for choosing Branson!

Introduction

This manual is arranged into several structured chapters which will help you find the information you may need to know to safely handle, install, set up, program, operate, and/or maintain this product. Please refer to the [Table of Contents](#) and/or the [Index](#) of this manual to find the information you may be looking for. In the event you require additional assistance or information, please contact our Product Support department (see [7.4 Service Centers](#) for information on how to contact them) or your local Branson representative.



Table of Contents

Chapter 1: Safety

| | |
|--|---|
| 1.1 Safety Requirements and Warnings | 2 |
| 1.2 Safety Precautions | 4 |
| 1.3 Warranty. | 6 |

Chapter 2: Introduction

| | |
|-------------------------------------|---|
| 2.1 How Ultrasonics Works | 8 |
| 2.2 Ultrasonic Baths | 9 |

Chapter 3: Delivery and Handling

| | |
|-----------------------------------|----|
| 3.1 Unpacking Your Unit | 12 |
|-----------------------------------|----|

Chapter 4: Technical Specifications

| | |
|--|----|
| 4.1 Model Name Definition | 14 |
| 4.2 Equipment Specifications. | 15 |
| 4.3 Temperature | 21 |
| 4.4 Cleaning Solutions | 22 |
| 4.5 Solution Effect on Metals. | 24 |

Chapter 5: Installation and Setup

| | |
|-----------------------------------|----|
| 5.1 Installing Your Unit. | 26 |
|-----------------------------------|----|

Chapter 6: Operation

| | |
|--|----|
| 6.1 Operating Your Ultrasonic Bath | 28 |
| 6.2 M Series | 29 |
| 6.3 MH Series | 32 |
| 6.4 CPX Series. | 35 |
| 6.5 CPXH Series | 40 |
| 6.6 Cleaning Methods. | 55 |

Chapter 7: Maintenance

| | |
|---|----|
| 7.1 Optimizing Your Ultrasonic Bath | 58 |
| 7.2 Troubleshooting | 59 |
| 7.3 Glass Slide Test | 61 |
| 7.4 Service Centers | 62 |
| 7.5 Information for Users on Disposal of Equipment. | 65 |



List of Figures

Chapter 1: Safety**Chapter 2: Introduction**

Figure 2.1 Unit with Digital Control, plus Heat and Timer9

Chapter 3: Delivery and Handling**Chapter 4: Technical Specifications****Chapter 5: Installation and Setup****Chapter 6: Operation**

Figure 6.1 M Series Controls30

Figure 6.2 MH Series Controls33

Figure 6.3 CPX Series Controls36

Figure 6.4 CPXH Series Controls41

Figure 6.5 Draining of Units 1800 and 280052

Figure 6.6 Direct Cleaning Method55

Figure 6.7 Indirect Cleaning Method55

Figure 6.8 Non Cleaning Application56

Chapter 7: Maintenance



List of Tables

Chapter 1: Safety

Chapter 2: Introduction

| | | |
|-----------|--------------------------------------|---|
| Table 2.1 | Ultrasonic Baths Available | 9 |
|-----------|--------------------------------------|---|

Chapter 3: Delivery and Handling

Chapter 4: Technical Specifications

| | | |
|------------|---|----|
| Table 4.1 | Model Name Definition | 14 |
| Table 4.2 | Equipment Specifications for North America Models | 15 |
| Table 4.3 | Equipment Specifications for Europe Models | 16 |
| Table 4.4 | Equipment Specifications for Japan Models | 17 |
| Table 4.5 | Equipment Specifications for China Models | 18 |
| Table 4.6 | Fuse Table for North America and Japan Models. | 19 |
| Table 4.7 | Fuse Table for Europe and China Models | 20 |
| Table 4.8 | Temperature | 21 |
| Table 4.9 | Alkaline Solution Strength and Uses | 22 |
| Table 4.10 | Chemicals Harmful to Your Tank | 23 |
| Table 4.11 | Solution Effects on Metals | 24 |

Chapter 5: Installation and Setup

Chapter 6: Operation

| | | |
|------------|--|----|
| Table 6.1 | Before you Begin | 29 |
| Table 6.2 | M Series Explanation of Controls | 30 |
| Table 6.3 | For initial cleaning solution degassing. | 30 |
| Table 6.4 | Treating Samples | 31 |
| Table 6.5 | Before you Begin | 32 |
| Table 6.6 | MH Series Explanation of Controls | 33 |
| Table 6.7 | For initial cleaning solution degassing. | 33 |
| Table 6.8 | Treating Samples | 34 |
| Table 6.9 | Before you begin | 35 |
| Table 6.10 | CPX Series explanation of controls | 36 |
| Table 6.11 | CPX Series LCD description. | 37 |
| Table 6.12 | Degassing your CPX Series unit | 38 |
| Table 6.13 | Treating Samples | 39 |
| Table 6.14 | Before you begin. | 40 |
| Table 6.15 | CPXH Series explanation of controls. | 41 |
| Table 6.16 | LCD Description for CPXH Series. | 45 |
| Table 6.17 | Degassing | 47 |
| Table 6.18 | Ultrasonic Operating Modes | 47 |
| Table 6.19 | Treating samples in Timed Sonics Mode | 48 |
| Table 6.20 | Treating Samples in Constant Sonics Mode | 49 |
| Table 6.21 | Treating Samples in Auto Mode | 50 |
| Table 6.22 | CPXH temperature calibration | 51 |
| Table 6.23 | Draining your ultrasonic bath. | 52 |
| Table 6.24 | Measuring the Solution Temperature | 54 |
| Table 6.25 | Solution Usage. | 54 |

Chapter 7: Maintenance

| | | |
|-----------|--|----|
| Table 7.1 | Tanks | 58 |
| Table 7.2 | Troubleshooting | 59 |
| Table 7.3 | Authorized Service Centers (North America) | 62 |
| Table 7.4 | Technical Support (North America) | 62 |
| Table 7.5 | Authorized Repair Representatives | 62 |
| Table 7.6 | Authorized Service Center/Technical Support (Europe) | 63 |
| Table 7.7 | Authorized Service Center/Technical Support (Asia) | 63 |

Chapter 1: Safety





| | | |
|------------|---|----------|
| 1.1 | Safety Requirements and Warnings | 2 |
| 1.2 | Safety Precautions | 4 |
| 1.3 | Warranty | 6 |

1.1 Safety Requirements and Warnings

This chapter contains an explanation of the different Safety Notice symbols and icons found both in this manual and on the product itself and provides additional safety information for Industrial Cleaning. This chapter also describes how to contact Branson for assistance.

1.1.1 Symbols Found in this Manual

These symbols used throughout this manual warrant special attention.


| | |
|---|---|
| DANGER | Indicates an immediate danger |
|  | If these risks are not avoided, death or severe injury will be the result. |
| WARNING | Indicates a possible danger |
|  | If these risks are not avoided, death or severe injury might result. |
| CAUTION | Indicates a possible danger |
|  | If these risks are not avoided, slight or minor injury might result. |
| NOTICE | Indicates a possible damaging situation |
|  | If this situation is not avoided, the system or something in its vicinity might get damaged. Application types and other important or useful information are emphasized. |

PN indicates Part Number.

Part(s) indicates your work piece or component to be cleaned.

1.1.2 Flammable chemicals

The use of flammable or explosive chemistry or materials in or around this machine is strictly prohibited. This machine is not designed for use with flammable/explosive chemistry nor should it be used in the same area as other equipment that uses flammable/explosive chemistry or materials. This machine is not intrinsically safe and can provide a source of ignition (from heat or electrical arc) for flammable or explosive chemicals or vapors which could cause a fire or explosion and may result in severe injury or death.

| WARNING | |
|---|--|
|  | NEVER use flammable or explosive material in the Ultrasonic Bath. |

1.2 Safety Precautions

Before using your Ultrasonic Bath, please read and thoroughly understand these safety precautions. Failure to follow them may result in serious personal injury or property damage.

To avoid electrical shock

- Do unplug from power source before filling or emptying the tank
- Do plug the unit into an appropriate grounded power socket
- Do connect the unit to a power supply using a properly sized overcurrent protection device. See label on the back of unit for information on current rating
- Do keep the control panel and the area around the unit clean and dry—wipe up solution which spills over the tank brim. Water and high voltage can cause electrical shock
- Do not operate the unit without proper grounding
- Do not remove the grounding prong on the line cord plug
- Do not disassemble your unit—high voltage inside the unit is dangerous
- Do not immerse the unit in water

To prevent personal and/or property damage

- Do use water-based solutions
- Do not ever use alcohol, gasoline or flammable solutions. Doing so could cause a fire or explosion and will void your warranty. Use only water-based solutions
- Do not ever use mineral acids. These could damage the tank
- Do not touch the stainless steel tank or cleaning solution—they may be hot
- Do not allow fluid temperature to exceed 70° C (158° F)
- Do not place your fingers or hands into the tank while the unit is operating. Doing so may cause discomfort and possible skin irritation. Avoid contact with solutions and provide adequate ventilation
- Do not use solutions containing chlorine bleach

To prevent damage to the unit

- Do change your solution regularly
- Do not cover vents on the cover
- Do not operate the unit dry
- Do not place parts or containers directly on the bottom of the tank; use a tray or wire to suspend items. Failure to comply may cause transducer damage and will void your warranty
- Do not allow the solution to drop more than 3/8 inch (1 cm) below the operating level line with heat or ultrasonics on. Failure to comply may cause transducer and/or heater damage and will void your warranty
- Turn off AC and heater switch before plugging/unplugging the line cord

Sound level and energy savings

- Do not operate the unit without a cover when possible
- The sound pressure released by the unit depends on the size of the bath and the application, but is less than 80 dBA when used with a cover
- To reduce the sound pressure it is recommended to use a cover while ultrasonics are activated and to switch the ultrasonics on with the bath loaded when possible

Insulation resistance test


Branson has taken all applicable measures to assure that manufactured units comply with insulation resistance requirements, as outlined by IEC 61010-1:2010 (Third edition). As per the Portable Appliance Testing (PAT) requirements, testing should be carried out by the user.

Before doing the insulation testing, however, please read the following pertinent information:


All of Branson's 220 V "C" and 230 V "E" units are equipped with Metal Oxide Varistors (MOVs) as the primary components to absorb overvoltages in the power line. MOVs disperse the over-voltage to the ground line and thus away from the equipment. Per their mode of operation, these MOVs will cause the insulation resistance test to fail when it is carried out at 500 V DC. As such, Branson recommends performing the test at a reduced voltage (250 V DC), as allowed by the Code of Practice for In-service Inspection and Testing of Electrical Equipment published by the IET (The Institution of Engineering and Technology), as this will prevent the MOVs from triggering and failing the test.

1.3 Warranty

For warranty information please reference the warranty section of Terms and Conditions found at: www.emerson.com/branson-terms-conditions.

| WARNING | General Warning |
|---|---|
|  | <ul style="list-style-type: none"> • Do not place parts or containers directly on the bottom of the tank; use a tray or wire to suspend items. Direct placement can cause the units to fail • Do not allow the solution to drop more than 3/8 inch (1 cm) below the operating level line • Do not ever use alcohol, gasoline or flammable solutions. Doing so could cause a fire or explosion. Use only water-based solutions • Do not use mineral acids. These could damage the tank |

First time cleaning - First experiment with one piece, then proceed with the remainder.

| CAUTION | General Warning |
|--|--|
|  | <p>Never clean novelty or inexpensive jewelry in the ultrasonic bath. The combination of heat and vibration may loosen a cement-held setting.</p> <p>Never clean gemstones such as emerald, amethyst, pearl, opal, coral, turquoise, peridot or lapis lazuli in the ultrasonic bath.</p> |

Solution level - Be sure to maintain solution level within 1/2 inch (1.3 cm) of the tank's "operating level" line. Surface activity can vary with liquid level.

Load size - It is faster and more efficient to run several small loads rather than a few big loads.

Placing items - Never allow items to sit on the bottom of the tank. Always place them in a tray or beaker or suspend in the solution.

Rinsing items - After cleaning, use a clean water bath to rinse away chemicals adhering to items.

Lubricating items - When necessary, re-lubricate items immediately after cleaning.

Drying items - Air drying at room temperature works for some items. Place parts requiring faster drying time under hot air blowers or in ovens.

Please call your local distributor if you have application questions.


Chapter 2: Introduction

| | | |
|------------|------------------------------------|----------|
| 2.1 | How Ultrasonics Works | 8 |
| 2.2 | Ultrasonic Baths | 9 |


2.1 How Ultrasonics Works

Ultrasonic sound is sound transmitted at frequencies generally beyond the range of human hearing. In your ultrasonic bath, ultrasonic sound (sonics) can be used for cleaning materials and parts, and for dissolving, homogenizing and degassing liquids. This is how it works:


- As the sound waves from the transducer radiate through the solution in the tank, they cause alternating high and low pressures in the solution



- During the low pressure stage, millions of microscopic bubbles form and grow. This process is called CAVITATION, meaning "formation of cavities"



- During the high pressure stage, the bubbles collapse, or "implode" releasing enormous amounts of energy
- For ultrasonic cleaning applications, these implosions act like an army of tiny scrub brushes. They work in all directions, attacking every surface and invading all recesses and openings



- This same energy can be used for other applications, such as liquid dissolving, homogenizations, and degassing

2.2 Ultrasonic Baths

This line of ultrasonic baths include five sizes:


Table 2.1 Ultrasonic Baths Available

| Model Number | Tank Capacity |
|--------------|----------------------|
| 1800 | 1/2 gal. (1.91 l) |
| 2800 | 3/4 gal. (2.81 l) |
| 3800 | 1-1/2 gal. (5.71 l) |
| 5800 | 2-1/2 gal. (9.51 l) |
| 8800 | 5-1/2 gal. (20.81 l) |

Each model is constructed using durable industrial style 40 kHz transducers. These provide increased ultrasonic power along with built in sweep frequency to ensure uniform ultrasonic activity throughout the bath. Models 1800 and 2800 have a molded dip in the left side of their rims to facilitate emptying of solution from the tank. Models 3800, 5800 and 8800 have built in drains and are supplied with tank drain kits. Each model can be purchased in four different configurations:

- With a Mechanical Timer (M)
- With a Mechanical Timer plus Heat (MH)
- With Digital Control and Timer (CPX)
- With Digital Control, plus Heat and Timer (CPXH)


Figure 2.1 Unit with Digital Control, plus Heat and Timer



When you first fill your unit, or refill it with fresh solution, use warm water for the solution. Turn on the heater (if available), turn on the ultrasonics (press the Sonics key or rotate the Timer), add the cover and the solution will heat quickly to temperature.

2.2.1 Accessories For Your Unit

As parts cannot be placed on the tank bottom, accessories include beaker positioning covers, solid and perforated insert trays, mesh baskets, beakers, and support racks.

| NOTICE | |
|---|---|
|  | Tank covers are included with every unit. |

Chapter 3: Delivery and Handling

3.1 Unpacking Your Unit..... 12

3.1 Unpacking Your Unit

Please check your unit and its carton carefully for any external or internal damage. If you find damage, contact your shipping carrier immediately, before contacting your distributor. Please retain your packaging for future use.

Chapter 4: Technical Specifications

| | | |
|------------|--|-----------|
| 4.1 | Model Name Definition | 14 |
| 4.2 | Equipment Specifications | 15 |
| 4.3 | Temperature | 21 |
| 4.4 | Cleaning Solutions | 22 |
| 4.5 | Solution Effect on Metals | 24 |

4.1 Model Name Definition

The name of the models determine the specifications of each unit. For example the CPX1800H-E:




Table 4.1 Model Name Definition

| Item | Stands For | Availability |
|------|----------------|---------------------------------------|
| 1 | Model | CPX: Digital |
| | | M: Mechanical |
| 2 | Tank Capacity | 1800: 1/2 gal (1.91 l) |
| | | 2800: 3/4 gal (2.81 l) |
| | | 3800: 1-1/2 gal (5.71 l) |
| | | 5800: 2-1/2 gal (9.51 l) |
| 3 | Heater | Blank: No heater |
| | | H: Heater |
| 4 | Region/voltage | Blank: North America (120 VAC) |
| | | E: Europe (230 VAC) |
| | | J: Japan (100 VAC) |
| | | C: China (220 VAC) |

- All models have a frequency of 40 kHz
- In CPXH models, the temperature readout accuracy is $\pm 3^{\circ}\text{C}$ ($\pm 5.4^{\circ}\text{F}$)
- Models available for 120 V $\pm 10\%$, 50/60 Hz and 220 V $\pm 10\%$, 50/60 Hz operation
- All 120 V units have CSA/UL or equivalent approval and comply with FCC regulations
- All 220-230 V units meet CE standards
- All units have a ground leakage current less than .50 ma
- Operating ambient temperature is from 5° C to 40° C (41° F to 104° F)

4.2 Equipment Specifications

Table 4.2 Equipment Specifications for North America Models

| Model Name | Tank Capacity | Tank Size (inches) | Overall Size (inches) | Weight | Max Sonics Power | Heater Power | Max. Draw Power Req. |
|------------|------------------------|-------------------------------|-------------------------------|--------------------|------------------|--------------|----------------------|
| M1800 | 1/2 gal. (1.91 l) | L: 9.9 W: 5.5 H: 4 | L: 9.9 W: 12 H: 11.9 | 9 lb (4 kg) | 70 | 0 | 90 |
| M1800H | | | | | 70 | 60 | 150 |
| CPX1800H | | | | | 70 | 60 | 150 |
| CPX1800 | | | | | 70 | 0 | 90 |
| M2800 | 3/4 gal (2.81 l) | L: 9.5 W: 5.5 H: 4 | L: 13.3 W: 12 H: 11.9 | 10 lb (4.5 kg) | 110 | 0 | 130 |
| M2800H | | | | | 110 | 90 | 250 |
| CPX2800H | | | | | 110 | 90 | 250 |
| CPX2800 | | | | | 110 | 0 | 250 |
| M3800 | 1-1/2 gal (5.71 l) | L: 11.5 W: 6 H: 6 | L: 15.6 W: 12.5 H: 14.8 | 14 lb (6.4 kg) | 110 | 0 | 130 |
| M3800H | | | | | 110 | 180 | 350 |
| CPX3800H | | | | | 110 | 180 | 350 |
| CPX3800 | | | | | 110 | 0 | 130 |
| M5800 | 2-1/2 gal (9.51 l) | L: 11.5 W: 9.5 H: 6 | L: 15.6 W: 15.8 H: 14.9 | 16 lb (7.3 kg) | 160 | 0 | 180 |
| M5800H | | | | | 160 | 280 | 490 |
| CPX5800H | | | | | 160 | 280 | 490 |
| CPX5800 | | | | | 160 | 0 | 180 |
| M8800 | 5-1/2 gal (20.81 l) | L: 19.5 W: 11.5 H: 15.4 | L: 23.5 W: 18.3 H: 15.4 | 28 lb (12.7 kg) | 280 | 0 | 320 |
| M8800H | | | | | 280 | 560 | 930 |
| CPX8800H | | | | | 280 | 560 | 930 |
| CPX8800 | | | | | 280 | 0 | 320 |

Table 4.3 Equipment Specifications for Europe Models

| Model Name | Tank Capacity | Tank Size (mm) | Overall Size (mm) | Weight | Max Sonics Power | Heater Power | Max. Draw Power Req. |
|------------|------------------------|----------------------------|----------------------------|--------------------|------------------|--------------|----------------------|
| M1800-E | 1.91 l (1/2 gal) | L: 150 W: 140 H: 100 | L: 251 W: 305 H: 302 | 5.4 kg (12 lb) | 70 | 0 | 90 |
| M1800H-E | | | | | 70 | 60 | 150 |
| CPX1800H-E | | | | | 70 | 60 | 150 |
| CPX1800-E | | | | | 70 | 0 | 90 |
| M2800-E | 2.81 l (3/4 gal) | L: 240 W: 140 H: 100 | L: 338 W: 305 H: 302 | 6.8 kg (15 lb) | 110 | 0 | 130 |
| M2800H-E | | | | | 110 | 110 | 250 |
| CPX2800H-E | | | | | 110 | 110 | 250 |
| CPX2800-E | | | | | 110 | 0 | 250 |
| M3800-E | 5.71 l (1-1/2 gal) | L: 290 W: 150 H: 150 | L: 396 W: 318 H: 302 | 8.2 kg (18 lb) | 110 | 0 | 130 |
| M3800H-E | | | | | 110 | 215 | 350 |
| CPX3800H-E | | | | | 110 | 215 | 350 |
| CPX3800-E | | | | | 110 | 0 | 350 |
| M5800-E | 9.51 l (2-1/2 gal) | L: 290 W: 240 H: 150 | L: 396 W: 401 H: 378 | 9.5 kg (21 lb) | 160 | 0 | 180 |
| M5800H-E | | | | | 160 | 300 | 490 |
| CPX5800H-E | | | | | 160 | 300 | 490 |
| CPX5800-E | | | | | 160 | 0 | 180 |
| M8800-E | 20.81 l (5-1/2 gal) | L: 495 W: 290 H: 150 | L: 597 W: 465 H: 391 | 16.3 kg (36 lb) | 280 | 0 | 320 |
| M8800H-E | | | | | 280 | 600 | 930 |
| CPX8800H-E | | | | | 280 | 600 | 930 |
| CPX8800-E | | | | | 280 | 0 | 320 |

Table 4.4 Equipment Specifications for Japan Models

| Model Name | Tank Capacity | Tank Size (mm) | Overall Size (mm) | Weight | Max Sonics Power | Heater Power | Max. Draw Power Req. |
|------------|------------------------|----------------------------|----------------------------|--------------------|------------------|--------------|----------------------|
| M1800-J | 1.91 l (1/2 gal) | L: 150 W: 140 H: 100 | L: 251 W: 305 H: 302 | 4 kg (9 lb) | 70 | 0 | 90 |
| M1800H-J | | | | | 70 | 45 | 135 |
| CPX1800H-J | | | | | 70 | 45 | 140 |
| CPX1800-J | | | | | 70 | 0 | 90 |
| M2800-J | 2.81 l (3/4 gal) | L: 240 W: 140 H: 100 | L: 338 W: 305 H: 302 | 4.5 kg (10 lb) | 110 | 0 | 130 |
| M2800H-J | | | | | 110 | 65 | 205 |
| CPX2800H-J | | | | | 110 | 65 | 205 |
| CPX2800-J | | | | | 110 | 0 | 130 |
| M3800-J | 5.71 l (1-1/2 gal) | L: 290 W: 150 H: 150 | L: 396 W: 318 H: 376 | 6.4 kg (14 lb) | 110 | 0 | 130 |
| M3800H-J | | | | | 110 | 130 | 275 |
| CPX3800H-J | | | | | 110 | 130 | 280 |
| CPX3800-J | | | | | 110 | 0 | 130 |
| M5800-J | 9.5 l (2-1/2 gal) | L: 290 W: 240 H: 150 | L: 396 W: 401 H: 378 | 7.3 kg (16 lb) | 160 | 0 | 180 |
| M5800H-J | | | | | 160 | 200 | 405 |
| CPX5800H-J | | | | | 160 | 200 | 410 |
| CPX5800-J | | | | | 160 | 0 | 180 |
| M8800-J | 20.81 l (5-1/2 gal) | L: 495 W: 290 H: 150 | L: 597 W: 465 H: 391 | 12.7 kg (28 lb) | 280 | 0 | 320 |
| M8800H-J | | | | | 280 | 400 | 755 |
| CPX8800H-J | | | | | 280 | 400 | 760 |
| CPX8800-J | | | | | 280 | 0 | 320 |

Table 4.5 Equipment Specifications for China Models

| Model Name | Tank Capacity | Tank Size (mm) | Overall Size (mm) | Weight | Max Sonics Power | Heater Power | Max. Draw Power Req. |
|------------|------------------------|----------------------------|----------------------------|--------------------|------------------|--------------|----------------------|
| M1800-C | 1.91 l (1/2 gal) | L: 150 W: 140 H: 100 | L: 251 W: 305 H: 302 | 5.4 kg (12 lb) | 70 | 0 | 90 |
| M1800H-C | | | | | 70 | 55 | 145 |
| CPX1800H-C | | | | | 70 | 55 | 145 |
| CPX1800-C | | | | | 70 | 0 | 90 |
| M2800-C | 2.81 l (3/4 gal) | L: 240 W: 140 H: 100 | L: 338 W: 305 H: 302 | 6.8 kg (15 lb) | 110 | 0 | 130 |
| M2800H-C | | | | | 110 | 105 | 250 |
| CPX2800H-C | | | | | 110 | 105 | 250 |
| CPX2800-C | | | | | 110 | 0 | 130 |
| M3800-C | 5.71 l (1-1/2 gal) | L: 290 W: 150 H: 150 | L: 396 W: 318 H: 376 | 8.2 kg (18 lb) | 110 | 0 | 130 |
| M3800H-C | | | | | 110 | 205 | 350 |
| CPX3800H-C | | | | | 110 | 205 | 350 |
| CPX3800-C | | | | | 110 | 0 | 130 |
| M5800-C | 9.5 l (2-1/2 gal) | L: 290 W: 240 H: 150 | L: 396 W: 401 H: 378 | 9.5 kg (21 lb) | 160 | 0 | 180 |
| M5800H-C | | | | | 160 | 285 | 490 |
| CPX5800H-C | | | | | 160 | 285 | 490 |
| CPX5800-C | | | | | 160 | 0 | 180 |
| M8800-C | 20.81 l (5-1/2 gal) | L: 495 W: 290 H: 150 | L: 597 W: 465 H: 391 | 16.3 kg (36 lb) | 280 | 0 | 320 |
| M8800H-C | | | | | 280 | 560 | 930 |
| CPX8800H-C | | | | | 280 | 560 | 930 |
| CPX8800-C | | | | | 280 | 0 | 320 |

Table 4.6 Fuse Table for North America and Japan Models

| Model Name | Fuse 1 | Fuse 2 | Fuse 3 |
|-----------------------|-------------|--------------|-----------|
| M1800 / M1800-J | 250 V, 2A | 250 V, 1.6A | 250 V, 1A |
| M1800H / M1800H-J | | | |
| CPX1800H / CPX1800H-J | | | |
| CPX1800 / CPX1800-J | 250 V, 1.6A | | |
| M2800 / M2800-J | 250 V, 2.5A | | |
| M2800H / M2800H-J | | | |
| CPX2800H / CPX2800H-J | | | |
| CPX2800 / CPX2800-J | 250 V, 1.6A | | |
| M3800 / M3800-J | 250 V, 2.5A | | |
| M3800H / M3800H-J | | | |
| CPX3800H / CPX3800H-J | | | |
| CPX3800 / CPX3800-J | 250 V, 1.6A | | |
| M5800 / M5800-J | 250 V, 5A | 250 V, 2A | |
| M5800H / M5800H-J | | | |
| CPX5800H / CPX5800H-J | | | |
| CPX5800 / CPX5800-J | 250 V, 2.5A | 250 V, 2.5A | |
| M8800 / M8800-J | 250 V, 10A | 250 V, 3.15A | |
| M8800H / M8800H-J | | | |
| CPX8800H / CPX8800H-J | | | |
| CPX8800 / CPX8800-J | 250 V, 5A | | |

Table 4.7 Fuse Table for Europe and China Models


| Model Name | Fuse 1 | Fuse 2 | Fuse 3 |
|-------------------------|-------------|-------------|-------------|
| M1800-E / M1800-C | 250 V, 1.6A | 250 V, 1.6A | 250 V, 1A |
| M1800H-E / M1800H-C | | | |
| CPX1800H-E / CPX1800H-C | | | |
| CPX1800-E / CPX1800-C | | | |
| M2800-E / M2800-C | 250 V, 2.5A | | |
| M2800-E / M2800H-C | | | |
| CPX2800H-E / CPX2800H-C | | | |
| CPX2800-E / CPX2800-C | 250 V, 1.6A | | |
| M3800-E / M3800-C | 250 V, 2.5A | | |
| M3800H-E / M3800H-C | | | |
| CPX3800H / CPX3800H-C | | | |
| CPX3800-E / CPX3800-C | 250 V, 1.6A | | |
| M5800-E / M5800-C | 250 V, 5A | 250 V, 2A | |
| M5800H-E / M5800H-C | | | |
| CPX5800H-E / CPX5800H-C | | | |
| CPX5800-E / CPX5800-C | 250 V, 1.6A | 250 V, 1.6A | |
| M8800-E / M8800-C | 250 V, 5A | 250 V, 2A | |
| M8800H-E / M8800H-C | | | |
| CPX8800H-E / CPX8800H-C | | | |
| CPX8800-E / CPX8800-C | | | 250 V, 2.5A |

4.3 Temperature

Table 4.8 Temperature

| Item | Function |
|---|---|
| Heater | The heater may cause some discoloration of the tank. This is normal and will not affect the performance of the unit. |
| Solution | The fastest method to heat your ultrasonic bath is to fill with warm solution, use heat, ultrasonics (which also adds heat), and a cover. |
| Over Temperature Protection (CPXH only) | If Max Temperature of 75° C is reached, Ultrasonics will pause for a period of time until the temperature in the bath falls back down to 69° C (the max set point). Once the temperature falls back down to 69° C, Ultrasonics will resume. This will apply for both continuous/infinity mode, as well as timed mode. In timed mode, the timer will pause while Ultrasonics is off, and will resume once Ultrasonics resumes. |

4.4 Cleaning Solutions

| CAUTION | General Warning |
|---|---|
|  | <p>Do not use alcohol, gasoline, bleach, mineral acids, solutions with a flash point, semi-aqueous or combustible liquids in ultrasonic tanks, or you will void the warranty. Only use non-flammable solutions and water-based solutions.</p> |

4.4.1 Solution Types

Water-based solutions can be either slightly acidic or alkaline. They include detergents, soaps and industrial cleaners designed to remove specific soils.

Acidic water-based solutions: remove rust, tarnish or scale. They range from mild solutions that remove tarnish, to concentrated, inhibited acidic solutions that remove investment plaster, milk-stone, zinc oxide and rust from steel and cast iron as well as smut and heat-treat scale from hardened steel.

Alkaline water-based solutions: include carbonates, silicates and caustics. These cause emulsifying action, which keeps soil from redepositing on the cleaned surface, and improves cleaning action in hard water.

Table 4.9 Alkaline Solution Strength and Uses

| Alkaline Strength | Removes |
|-------------------|---|
| Mild | Light oils and greases, cutting oils and coolant compounds. |
| Mild to Strong | Heavy grease and oils, waxes, vegetable oils, inks, wax or fat-base buffing and polishing compounds, milk residues and carbohydrates. |
| Heavy-duty | Mill scale, heat-treat scale, corrosion or oxides. |

Change the cleaning solution periodically. Cleaning solutions can become contaminated with suspended soil particles which coat the tank bottom. This coating dampens the ultrasonic action and reduces cleaning efficiency. Certain solutions will cavitate better than others. Contact your local distributor for further information.

Heat and cavitation increase the chemical activity of cleaning solutions. Some materials may be damaged by this stronger chemical action. When in doubt, test run samples of items to be cleaned.

Caustic solutions: used to remove rust from steels, metal alloy corrosion and a variety of tenacious soils.

4.4.2 Solution Amounts

Solution amounts may vary. The amount you use depends on the detergent and the type of soil to be removed. Follow instructions on the solution container and refer to the table below for the effects of solutions on metals.

4.4.3 Chemicals Harmful to Your Tank

The following chemicals will harm your ultrasonic tank and the action of ultrasonics and higher operating temperatures will increase their chemical activity. Do not use these or similar chemicals directly or in dilution in your ultrasonic tank or you will void your warranty.

Table 4.10 Chemicals Harmful to Your Tank

| Harmful Chemicals | | |
|----------------------|---------------------|-----------------------|
| Acetophenone | Chloracetic Acid | Hydrocyanic Acid |
| Aluminum Chloride | Chloric Acid | Hydrofluoric Acid |
| Aluminum Flouride | Chlorine, Anhydrous | Hydroflousilicic Acid |
| Aluminum Sulphate | Chromic Acid | Iodoform |
| Ammonium Bifluoride | Copper Chloride | Mercuric Chloride |
| Ammonium Chloride | Copper Fluoborate | Muriatic Acid |
| Ammonium Hydroxide | Ethyl Chloride | Phosphoric (crude) |
| Amyl Chloride | Ferric Chloride | Sodium Hypochlorite |
| Antimony Trichloride | Ferrous Chloride | Potassium Chloride |
| Aqua Regia | Ferris Sulfate | Stannic Chloride |
| Bromine | Fluoboric Acid | Stannous Chloride |
| Calcium Bisulfate | Fluorine | Sulfur Chloride |
| Calcium Bisulfite | Hydrobromic Acid | Sulfuric Acid |
| Calcium Hypochloride | Hydrochloric Acid | Zinc Chloride |

4.5 Solution Effect on Metals


Table 4.11 Solution Effects on Metals

| Cleaning Agent | Steel | Brass | Aluminum | Magnesium | Zinc | S. Steel Copper | Tin |
|--------------------------|-------------|--------------|---------------|------------|---------|-----------------|---------------|
| Optical (1) | none | none | none | none** | none** | none | none** |
| Jewelry (1) | none | none | none | none | none | none | none |
| Buffing (1) Compound | none | slight stain | none | none | attacks | none | none |
| Oxide (2) remover | slight etch | none | slight attack | attacks | attacks | none | none |
| Electronic cleaner (1) | none | none | slight attack | none | none | none | none |
| General (1) purpose | none | none | slight attack | none | none | none | none |
| Industrial strength (1) | none | none | slight attack | none | none | none | none |
| Metal (1) cleaner 1 | none | none | none | none | none | none | none |
| Metal (1) cleaner 2 | none | none | slight attack | none | none | none | none |
| Metal (1) cleaner 3 | none | none | none | none | none | none | none |
| Liquid rust (3) stripper | none | none | attacks*** | attacks*** | attacks | none | slight attack |
| GP (1) Powder | none | none | none | none | none | none | none |

* Contact distributor for cleaning agent availability outside the US.

** No effect if solution temperature is less than 60° C (140° F).

(1) = Alkaline; (2) = Acidic; and (3) = Caustic.

| WARNING | General Warning |
|---|---|
|  | <p>*** Free hydrogen may be released if solution comes in contact with reactive metals.</p> |

Chapter 5: Installation and Setup

5.1 Installing Your Unit 26

5.1 Installing Your Unit

Check the plate on the back of the unit for correct power requirements. Position your unit within easy reach of a standard grounded electrical outlet. Do not place the unit on a circuit which could become overloaded. If your unit does not operate correctly, first refer to [7.2 Troubleshooting](#) for possible causes, or contact an authorized service center listed at the end of this manual for additional information.

Chapter 6: Operation

| | | |
|------------|---|-----------|
| 6.1 | Operating Your Ultrasonic Bath | 28 |
| 6.2 | M Series | 29 |
| 6.3 | MH Series | 32 |
| 6.4 | CPX Series | 35 |
| 6.5 | CPXH Series | 40 |
| 6.6 | Cleaning Methods | 55 |

6.1 Operating Your Ultrasonic Bath

If this is the first time you are using the ultrasonic bath, please read this whole section before operating your unit.

6.2 M Series



| CAUTION | General Warning |
|---|---|
|  | <ul style="list-style-type: none"> Do not place parts or containers directly on the bottom of the tanks; use a tray or wire to suspend items. Direct placement can cause the units to fail Do not allow the solution to drop more than 3/8 inch (1 cm) below the operating level line with ultrasonics on Do not ever use alcohol, gasoline or flammable solutions. Doing so could cause a fire or explosion. Use only water-based solutions Do not ever use mineral acids. These could damage the tank <p>Failure to comply with these cautions will void your warranty.</p> |

Table 6.1 Before you Begin

| Step | Action |
|------|--|
| 1 | Select your cleaning solution (refer to 4.5 Solution Effect on Metals). |
| 2 | Allowing for the volume of the parts you will be cleaning and leaving room for cleaning solution, fill the tank with warm tap water to the operating level line. |
| 3 | Add cleaning solution to the tank water. |
| 4 | Plug the unit into a grounded outlet. |
| 5 | For maximum efficiency, refer to 7.1 Optimizing Your Ultrasonic Bath , before proceeding. |

| NOTICE | |
|---|--|
|  | <p>If this is the first time you are running the unit, or if you have changed cleaning solution, you must degas the solution. If not, skip to 6.2.3 Cleaning Items (Treating Samples).</p> |

6.2.1 M Series Explanation of Controls

Figure 6.1 M Series Controls



Table 6.2 M Series Explanation of Controls

| Control | Function |
|-------------------|---|
| Main Power Switch | <p>This switch is located on the back of the unit, next to the power cord receptacle.</p> <ul style="list-style-type: none"> Press the on side to power on the unit Press the off side to power off the unit <p>When operating the unit, normally leave the Main Power switch in the on position, and use the Timer Knob to activate ultrasonics.</p> |
| Timer Knob | <p>Activates ultrasonics and sets time.</p> <ul style="list-style-type: none"> Turn clockwise for timed operation (0–60 minutes) Turn counterclockwise to the HOLD position for continuous operation Turn to the zero position to turn unit Off |

6.2.2 Degassing

For initial cleaning solution degassing.

Table 6.3 For initial cleaning solution degassing.

| Step | Action |
|------|---|
| 1 | Turn Main Power switch on. |
| 2 | <p>Turn the Timer Knob clockwise to 5–10 and let the unit run to allow the solution to “degas”.</p> <p>NOTICE Refer to 6.5.11 Solution for information on degassing.</p> |

6.2.3 Cleaning Items (Treating Samples)


| NOTICE | |
|---|---|
|  | <p>To stop ultrasonics at any time, turn the Timer Knob to the zero position.</p> |

Table 6.4 Treating Samples

| Step | Action |
|------|--|
| 1 | Turn Main Power switch on. |
| 2 | Turn the Timer Knob clockwise to set the amount of time (0 - 60 minutes) you wish the items to be cleaned. Turn the Timer Knob counterclockwise to the Hold position for continuous operation. |
| 3 | Place the items into a basket, perforated tray, or beakers in a positioning cover. |
| 4 | If using beakers or a solid tray, add cleaning solution to beakers or tray to cover the items. |
| 5 | Slowly lower the tray or beakers into the tank. Do not allow items to contact the tank bottom. |
| 6 | When items are clean, slowly remove them from the tank. |
| 7 | Rinse the clean items with clean water and dry them, if necessary. |

6.3 MH Series



| CAUTION | General Warning |
|---|---|
|  | <ul style="list-style-type: none"> Do not place parts or containers directly on the bottom of the tanks; use a tray or wire to suspend items. Direct placement can cause the units to fail Do not allow the solution to drop more than 3/8 inch (1 cm) below the operating level line with heat or ultrasonics on Do not ever use alcohol, gasoline or flammable solutions. Doing so could cause a fire or explosion. Use only water-based solutions Do not ever use mineral acids. These could damage the tank <p>Failure to comply with these cautions will void your warranty.</p> |

Table 6.5 Before you Begin

| Step | Action |
|------|--|
| 1 | Select your cleaning solution (refer to 4.5 Solution Effect on Metals). |
| 2 | Allowing for the volume of the parts you will be cleaning and leaving room for cleaning solution, fill the tank with warm tap water to the operating level line. |
| 3 | Add cleaning solution to the tank water. |
| 4 | Plug the unit into a grounded outlet. |
| 5 | For maximum efficiency, refer to 7.1 Optimizing Your Ultrasonic Bath , before proceeding. |

| NOTICE | |
|---|--|
|  | <p>If this is the first time you are running the unit, or if you have changed cleaning solution, you must degas the solution. If not, skip to 6.3.3 Cleaning Items (Treating Samples).</p> |

6.3.1 MH Series Explanation of Controls

Figure 6.2 MH Series Controls



Table 6.6 MH Series Explanation of Controls

| Control | Function |
|-------------------|---|
| Main Power Switch | <p>This switch is located on the back of the unit, next to the power cord receptacle.</p> <ul style="list-style-type: none"> Press the on side to power on the unit Press the off side to power off the unit <p>When operating the unit, normally leave the Main Power switch in the (on) position, and use the Timer Knob to activate ultrasonics.</p> |
| Heat Switch | <p>Activates heat to 60° C (140° F) maximum.</p> <p>NOTICE Refer to 4.3 Temperature for further information on temperature.</p> |
| Timer Knob | <p>Activates ultrasonics and sets time.</p> <ul style="list-style-type: none"> Turn clockwise for timed operation (0–60 minutes) Turn counterclockwise to the HOLD position for continuous operation <p>Turn to the zero position to turn unit Off.</p> |

6.3.2 Degassing

Table 6.7 For initial cleaning solution degassing.

| Step | Action |
|------|---|
| 1 | Turn Main Power switch on. |
| 2 | Turn Heat switch on. |
| 3 | <p>Turn the Timer Knob clockwise to 5–10 and let the unit run to allow the solution to “degas”.</p> <p>NOTICE Refer to 6.5.11 Solution for information on degassing.</p> |

6.3.3 Cleaning Items (Treating Samples)


| NOTICE | |
|---|---|
|  | <p>To stop ultrasonics at any time, turn the Timer Knob to the zero position.</p> |

Table 6.8 Treating Samples

| Step | Action |
|------|--|
| 1 | Turn Main Power switch on. |
| 2 | Turn the Timer Knob clockwise to set the amount of time (0–60 minutes) you wish the items to be cleaned. Turn the Timer Knob counterclockwise to the HOLD position for continuous operation. |
| 3 | Place the items into a basket, perforated tray, or beakers in a positioning cover. |
| 4 | If using beakers or a solid tray, add cleaning solution to beakers or tray to cover the items. |
| 5 | Slowly lower the tray or beakers into the tank. Do not allow items to contact the tank bottom. |
| 6 | When items are clean, slowly remove them from the tank. |
| 7 | Rinse the clean items with clean water and dry them, if necessary. |

6.4 CPX Series



| CAUTION | General Warning |
|---|---|
|  | <ul style="list-style-type: none"> Do not place parts or containers directly on the bottom of the tanks; use a tray or wire to suspend items. Direct placement can cause the units to fail Do not allow the solution to drop more than 3/8 inch (1 cm) below the operating level line with ultrasonics on Do not ever use alcohol, gasoline or flammable solutions. Doing so could cause a fire or explosion. Use only water-based solutions Do not ever use mineral acids. These could damage the tank <p>Failure to comply with these cautions will void your warranty.</p> |

Table 6.9 Before you begin

| Step | Action |
|------|--|
| 1 | Select your cleaning solution (refer to 4.5 Solution Effect on Metals). |
| 2 | Allowing for the volume of the parts you will be cleaning and leaving room for cleaning solution, fill the tank with warm tap water to the operating level line. |
| 3 | Add cleaning solution to the tank water. |
| 4 | Plug the unit into a grounded outlet. |
| 5 | For maximum efficiency, refer to 7.1 Optimizing Your Ultrasonic Bath , before proceeding. |

| NOTICE | |
|---|--|
|  | <p>If this is the first time you are running the unit, or if you have changed cleaning solution, you must degas the solution. If not, skip to 6.4.4 Cleaning Items (Treating Samples).</p> |

6.4.1 CPX Explanation of Controls

Figure 6.3 CPX Series Controls



Table 6.10 CPX Series explanation of controls





| Control | Function |
|---|---|
| Main Power Switch | <p>This switch is located on the back of the unit, next to the power cord receptacle.</p> <ul style="list-style-type: none"> • Press the on side to power on the unit • Press the off side to power off the unit <p>When operating the unit, normally leave the Main Power switch in the on position, and use the On/Standby key to switch between the operating state and standby state.</p> |
|  | <p>On/Standby</p> <p>When the Main Power switch on the rear panel is in the on position, press to power on/off the unit.</p> |
|  | <p>Up/Down Keys</p> <p>Press to increase/decrease ultrasonic or degassing cycle time (hold for quick increments/decrements).</p> <p>Time values are circular, pressing Up from 99 minutes takes you to Constant Sonics Mode (Constant Sonics icon and "- -" display on the LCD) and then to 1 minutes. Pressing Down from 1 minutes takes you to Constant Sonics Mode and then to 99 minutes.</p> <p>During power-up, use to select high or low ultrasonic power output.</p> |

Table 6.10 CPX Series explanation of controls

| Control | Function |
|---|--|
|  | <p>Sonics</p> <p>Press to activate ultrasonics. If running in Timed Mode, a timer will begin to count down and ultrasonics will stop at 0 minutes. In Constant Sonics Mode (Constant sonics icon and "- -" on the display), timer has no function.</p> <p>Press sonics key again to deactivate ultrasonics.</p> <p>If running in Timed Mode, press Up and Down keys to adjust the ultrasonic cycle time (adjustable from 1 to 99 minutes).</p> |
|  | <p>Degas</p> <p>Press to degas the solution or to run a degas application. A default timer of 5 minutes will begin to count down and degassing will stop at 0 minutes.</p> <p>Press Degas key again to stop degassing the solution.</p> <p>During a degas cycle, press Up and Down keys to adjust the degas cycle time (adjustable from 1 to 99 minutes).</p> <p>NOTICE Refer to 6.5.11 Solution for information on degassing.</p> |

6.4.2 CPX Series LCD Description

Table 6.11 CPX Series LCD description







| Reference | Function |
|---|--|
|  | <p>Power Level</p> <p>Displayed for 15 s only during power-up, shows the current ultrasonic output power selection.</p> <p>Press the Sonics or Degas key to go into normal operating mode.</p> <p>Press Up or Down keys to change between high (HI) and low (LO) power ultrasonics.</p> |
|  | <p>Sonics/Degas Timer</p> <p>Displays the duration of a timed ultrasonic or degas cycle.</p> <p>Press Up and Down keys to adjust ultrasonic or degassing cycle time (adjustable from 1 to 99 minutes).</p> <p>In Constant Sonics Mode, "- -" is displayed.</p> |
|  | <p>Constant Sonics</p> <p>Indicates the unit is operating in Constant Sonics Mode.</p> <p>In Constant Sonics Mode, sonics will remain on until the Sonics key is pressed or the unit is turned off.</p> |

Table 6.11 CPX Series LCD description

| Reference | Function |
|--|--|
|  | <p>Sonics On</p> <p>Indicates sonics are active.</p> <p>If running in Timed Mode, ultrasonics will remain on until the timer reaches 0 minutes.</p> <p>In Constant Sonics Mode, ultrasonics will remain on until the Sonics key is pressed or the unit is turned off.</p> |
|  | <p>Degas On</p> <p>Indicates the unit is in Degas Mode.</p> <p>In Degas Mode, degassing will continue until the timer reaches 0 minutes.</p> <p>NOTICE Refer to 6.5.11 Solution for information on degassing.</p> |
|  | <p>Alarm</p> <p>Alarm Bell icon flashes when the unit encounters an abnormal operating condition.</p> <p>NOTICE Refer to 7.2 Troubleshooting for information on troubleshooting.</p> |

6.4.3 Degassing

For initial cleaning solution degassing.


| NOTICE | |
|---|--|
|  | <p>To stop degassing at any time, press the Degas key.</p> |

Table 6.12 Degassing your CPX Series unit

| Step | Action |
|------|---|
| 1 | Turn Main Power switch on. |
| 2 | Press the On/Standby key to turn on the unit. |
| 3 | <p>Press Degas key once to start the degas process.</p> <p>Default degas time is 5 minutes.</p> <p>If necessary, use Up/Down keys to alter degas time during a degas cycle.</p> <p>NOTICE Refer to 6.5.11 Solution for information on degassing.</p> |

Table 6.12 Degassing your CPX Series unit

| Step | Action |
|------|---|
| 4 | After completing the degas time, you are ready to set operating parameters. |

6.4.4 Cleaning Items (Treating Samples)


| NOTICE | |
|---|--|
|  | To stop ultrasonics at any time, press the Sonics key. |

Table 6.13 Treating Samples

| Step | Action |
|------|---|
| 1 | Turn Main Power switch on. |
| 2 | Press the On/Standby key to turn on the unit. |
| 3 | Set the amount of time you wish the items to be cleaned, or select Constant Sonics Mode: <ul style="list-style-type: none"> • Use Up/Down keys to increase/decrease cycle time (hold for quick increments/decrements) • Pressing Up key from 99 minutes or Down key from 1 minutes takes you to Constant Sonics Mode (Constant Sonics icon and "- -" display on the LCD screen) |
| 4 | Press the Sonics key to activate ultrasonics. |
| 5 | Place the items into a basket, perforated tray, or beakers in a positioning cover. |
| 6 | If using beakers or a solid tray, add cleaning solutions to beakers or tray to cover the items. |
| 7 | Slowly lower the tray or beakers into the tank. Do not allow items to contact the tank bottom. |
| 8 | When items are clean, slowly remove them from the tank. |
| 9 | Rinse clean items with clean, warm water and dry, if necessary. |

6.5 CPXH Series



| CAUTION | General Warning |
|---|---|
|  | <ul style="list-style-type: none"> • Do not place parts or containers directly on the bottom of the tanks; use a tray or wire to suspend items. Direct placement can cause the units to fail • Do not allow the solution to drop more than 3/8 inch (1 cm) below the operating level line with heat or ultrasonics on • Do not ever use alcohol, gasoline or flammable solutions. Doing so could cause a fire or explosion. Use only water-based solutions • Do not ever use mineral acids. These could damage the tank <p>Failure to comply with these cautions will void your warranty.</p> |

Table 6.14 Before you begin

| Step | Action |
|------|--|
| 1 | Select your cleaning solution (refer to 4.5 Solution Effect on Metals). |
| 2 | Allowing for the volume of the parts you will be cleaning and leaving room for cleaning solution, fill the tank with warm tap water to the operating level line. |
| 3 | Add cleaning solution to the tank water. |
| 4 | Plug the unit into a grounded outlet. |
| 5 | For maximum efficiency, refer to 7.1 Optimizing Your Ultrasonic Bath , before proceeding. |

| NOTICE | |
|---|---|
|  | <p>If this is the first time you are running the unit, or if you have changed cleaning solution, you must degas the solution. If not, skip to 6.5.5 Cleaning Items (Treating Samples) in Timed Sonics Mode.</p> |

6.5.1 CPXH Series Explanation of Controls

Figure 6.4 CPXH Series Controls



Table 6.15 CPXH Series explanation of controls




| Control | Function |
|---|--|
| Main Power Switch | <p>This switch is located on the back of the unit, next to the power cord receptacle.</p> <ul style="list-style-type: none"> Press the on side to power on the unit Press the off side to power off the unit. When operating the unit, normally leave the Main Power switch in the on position, and use the On/Standby key to switch between the operating state and standby state |
|  | <p>On/Standby</p> <p>When the Main Power switch on the rear panel is in the on position, press to power on/off the unit.</p> |
|  | <p>Up/Down Keys</p> <p>Press to increase/decrease ultrasonic or degassing cycle time (hold for quick increments/decrements).</p> <p>Time values are circular, pressing Up key from 99 minutes takes you to 1 minutes. Pressing Down key from 1 minutes takes you to 99 minutes.</p> <p>If the Fn key was pressed, use the Up/Down keys to adjust function settings.</p> |
|  | <p>Heat</p> <p>Press to turn heater on/off.</p> <p>Heater will shut off when set temperature is attained.</p> |

Table 6.15 CPXH Series explanation of controls




| Control | Function |
|---|---|
|  | <p>Sonics</p> <p>Press to activate ultrasonics. If running in Timed Mode, a timer will begin to count down and ultrasonics will stop at 0 minutes. In Constant Sonics Mode (Constant Sonics icon and “- - -” on the display), timer has no function.</p> <p>Press Sonics key again to deactivate ultrasonics.</p> <p>If running in Timed Mode, press Up and Down keys to adjust the ultrasonic cycle time (adjustable from 1 to 99 minutes).</p> |
|  | <p>Degas</p> <p>Press to degas the solution or to run a degas application. The degas timer will begin to count down from its current setting and degassing will stop at 0 minutes.</p> <p>Press Degas key again to stop degassing the solution.</p> <p>During a degas cycle, press Up and Down keys to adjust degas time (adjustable from 1 to 99 minutes).</p> <p>NOTICE Refer to 6.5.11 Solution for information on degassing.</p> |
|  | <p>Auto</p> <p>Press to begin an auto cycle. In Auto Mode, the following actions are carried out automatically by the controller:</p> <ul style="list-style-type: none"> • Heater is turned on to bring bath to set temperature • When set temperature is reached, ultrasonics are activated. The unit will abort the auto cycle and flash the AUTO icon if set temperature is not reached within a 120-minute period • When ultrasonics timer reaches 0 minutes, the auto cycle is finished <p>If at any point during an auto cycle the degas key is pressed, a degas cycle will begin. If ultrasonics has already started, the ultrasonics timer will restart after the degas period.</p> |

Table 6.15 CPXH Series explanation of controls







| Control | Function |
|---|---|
|  | <p>Press the Fn key to access less-frequently used function.</p> <ul style="list-style-type: none"> • Press 1x to Set Target Temp • Press 2x to Set Constant Sonics • Press 3x to Set Power Level • Press 4x Sel. Temp Units • Press 5x to Set Degas Time • Press 6x Ready State (Set Sonics Time) <p>The appropriate icon will flash to indicate which option is selected.</p> <p>If no key is pressed after 15 seconds the unit will save any changes and return to the Ready state.</p> <p>Press the Fn key again after making any changes to scroll through the rest of the options and return to the Ready state.</p> <p>The following options are available:</p> <ul style="list-style-type: none"> • Set Temperature (Fn 1x): Press Up and Down keys to increase/decrease the bath set temperature |
| | <div style="display: flex; align-items: center; justify-content: center; margin-bottom: 10px;"> <div style="border: 1px solid black; border-radius: 15px; padding: 10px 20px; margin-right: 10px;">Fn</div> <div style="font-size: 24px; margin-right: 10px;">(1x)</div> <div style="border: 1px dashed black; border-radius: 15px; padding: 5px 10px; margin-left: 20px;">Indicates Blinking</div> </div> <div style="border: 1px solid black; padding: 5px; margin-bottom: 10px;">  </div> <ul style="list-style-type: none"> • Timed/Constant Sonics (Fn 2x): Press Up key to select Constant Sonics Mode (Constant Sonics icon and "-- --" on the display) Press down key to select Timed Sonics Mode <div style="display: flex; align-items: center; justify-content: center; margin-bottom: 10px;"> <div style="border: 1px solid black; border-radius: 15px; padding: 10px 20px; margin-right: 10px;">Fn</div> <div style="font-size: 24px; margin-right: 10px;">(2x)</div> <div style="border: 1px dashed black; border-radius: 15px; padding: 5px 10px; margin-left: 20px;">Indicates Blinking</div> </div> <div style="border: 1px solid black; padding: 5px;">  </div> |

Table 6.15 CPXH Series explanation of controls

| Control | Function |
|---|--|
| <div data-bbox="277 1087 380 1220" style="background-color: #0056b3; color: white; padding: 5px; width: 40px; margin: 0 auto;">Fn</div> | <ul style="list-style-type: none"> • Full/Low Power (Fn 3x): Press up key to select full power ultrasonic output. Press down key to select low power ultrasonic output <div data-bbox="565 369 1170 520" style="display: flex; align-items: center; justify-content: center;"> <div style="border: 1px solid black; border-radius: 15px; padding: 10px; margin-right: 20px;">Fn (3x)</div> <div style="border: 1px dashed black; border-radius: 15px; padding: 10px; margin-left: 20px;">Indicates Blinking</div> </div> |
| |  |
| | <ul style="list-style-type: none"> • Temperature Units (Fn 4x): Press key to select Fahrenheit (°F). Press Down key to select Celsius (°C) <div data-bbox="565 915 1170 1066" style="display: flex; align-items: center; justify-content: center;"> <div style="border: 1px solid black; border-radius: 15px; padding: 10px; margin-right: 20px;">Fn (4x)</div> <div style="border: 1px dashed black; border-radius: 15px; padding: 10px; margin-left: 20px;">Indicates Blinking</div> </div> |
| |  |
| | <ul style="list-style-type: none"> • Degas Time (Fn 5x): Press Up/Down keys to increase/decrease degas time (hold for quick increments/decrements). Time values are circular, pressing Up from 99 minutes takes you to 1 minute. Pressing Down from 1 minute takes you to 99 minutes |
| | <div data-bbox="565 1520 1170 1671" style="display: flex; align-items: center; justify-content: center;"> <div style="border: 1px solid black; border-radius: 15px; padding: 10px; margin-right: 20px;">Fn (5x)</div> <div style="border: 1px dashed black; border-radius: 15px; padding: 10px; margin-left: 20px;">Indicates Blinking</div> </div>  |

6.5.2 CPXH Series LCD Description

Table 6.16 LCD Description for CPXH Series












| Item | Function |
|---|--|
|  | <p>Sonics/Degas Set Time</p> <p>Displays the set time for a timed ultrasonic or degas cycle.</p> <ul style="list-style-type: none"> Press Up and Down keys to adjust ultrasonic or degassing cycle time (adjustable from 1 to 99 minutes). In Constant Sonics Mode, "-- --" is displayed |
|  | <p>Set Temperature</p> <p>Displays the target temperature.</p> <p>Temperature units are indicated by the °F (for Fahrenheit) or °C (for Celsius) right of the Current Temperature icon.</p> |
|  | <p>Current Temperature</p> <p>Displays the current tank temperature as measured by the unit.</p> <p>Temperature units are indicated right of the icon as either °F (For Fahrenheit) or °C (for Celsius). Units can be switched using the Fn key. See Fn key description on Table 6.15.</p> |
|  | <p>Sonics/Degas Timer</p> <p>Displays the remaining time of a running timed ultrasonic or degas cycle.</p> <p>Press up and down keys to adjust ultrasonic or degassing cycle time (adjustable from 1 to 99 minutes).</p> |
|  | <p>Auto</p> <p>Indicates the unit is in Auto Mode. In Auto Mode, the following actions are carried out automatically by the controller:</p> <ul style="list-style-type: none"> Heater is turned on to bring bath to set temperature When set temperature is reached, ultrasonics are activated. The unit will abort the auto cycle and flash the icon if set temperature is not reached within a 120-minute period When ultrasonics timer reaches 0 minutes, the auto cycle is finished <p>If at any point during an auto cycle the degas key is pressed, a degas cycle will begin. If ultrasonics has already started, the ultrasonics timer will restart after the degas period.</p> |
|  | <p>Constant Sonics</p> <p>Indicates the unit is operating in Constant Sonics Mode. In Constant Sonics Mode, ultrasonics will remain on until the Sonics key is pressed or the unit is turned off.</p> |

Table 6.16 LCD Description for CPXH Series

| Item | Function |
|---|--|
|  | <p>Heat</p> <p>Indicates the heater is on. Heater will shut off when set temperature is attained.</p> |
|  | <p>Sonics On</p> <p>Indicates sonics are active.</p> <p>If running in Timed Mode, ultrasonics will remain on until the timer reaches 0 minutes.</p> <p>In Constant Sonics Mode, ultrasonics will remain on until the Sonics key is pressed or the unit is turned off.</p> |
|  | <p>Degas On</p> <p>Indicates the unit is in Degas Mode.</p> <p>In Degas Mode, degassing will continue until the timer reaches 0 minutes.</p> <p>NOTICE</p> <p>Refer to 6.5.11 Solution for information on degassing.</p> |
|  | <p>Power Level</p> <p>Indicates the ultrasonic power output selection:</p> <ul style="list-style-type: none"> • Four bars indicate high power ultrasonics • Two bars indicate low power ultrasonics |
|  | <p>Alarm</p> <p>Alarm Bell icon flashes when the unit encounters an abnormal operating condition.</p> <p>NOTICE</p> <p>Refer to 7.2 Troubleshooting for information on troubleshooting.</p> |

6.5.3 Degassing

For initial cleaning solution degassing.


| NOTICE | |
|---|--|
|  | <p>To stop degassing at any time, press the Degas key.</p> |

Table 6.17 Degassing

| Step | Action |
|------|---|
| 1 | Turn Main Power switch on. |
| 2 | Press the On/Standby key to turn on the unit. |
| 3 | Default degas time is 5 minutes. To change the degas time, press the Fn key until the Degas icon appears and the Set Time icon flashes. Then press the Up/Down keys to change the degas time. |
| 4 | Press Degas key once to start the degas process. If necessary, use up/down keys to alter degas time during a degas cycle. NOTICE Refer to 6.5.11 Solution for information on degassing. |
| 5 | After completing the degas time, you are ready to set operating parameters. |

6.5.4 Ultrasonics Operating Modes

Table 6.18 Ultrasonic Operating Modes

| Mode | Action |
|-----------------|---|
| Timed Sonics | In Timed Sonics Mode, a timer will begin to count down and ultrasonics will remain on until the timer reaches 0 minutes. For instructions, see 6.5.5 Cleaning Items (Treating Samples) in Timed Sonics Mode . |
| Constant Sonics | In Constant Sonics Mode ultrasonics will remain on until the Sonics key is pressed or power is turned off to the unit. For instructions, see 6.5.6 Cleaning Items (Treating Samples) in Constant Sonics Mode . |
| Auto | In Auto mode ultrasonics will start once set temperature is attained. Ultrasonics will remain on until the timer reaches 0 minutes. For instructions, see 6.5.7 Cleaning Items (Treating Samples) in Auto Mode . |

6.5.5 Cleaning Items (Treating Samples) in Timed Sonics Mode


| NOTICE | |
|---|--|
|  | To stop ultrasonics at any time, press the Sonics key. |

Table 6.19 Treating samples in Timed Sonics Mode

| Step | Action |
|------|---|
| 1 | Turn Main Power on. |
| 2 | Press the On/Standby key to turn on the unit. |
| 3 | If necessary, degas the liquid. See 6.5.3 Degassing for instructions. |
| 4 | Set the amount of time you wish the items to be cleaned: Use up/down keys to increase/decrease cycle time (hold for quick increments/decrements). |
| 5 | Set the tank temperature: <ul style="list-style-type: none"> • Press the Fn key until the Set Temperature icon flashes • Press the Up/Down keys to alter the setting to the tank temperature you wish to reach • Press the Heat key once to activate heat. The Heat icon appears <div style="background-color: #003366; color: white; padding: 2px; margin-top: 5px;">NOTICE</div> Units can be switched between °F or °C using the Fn key. See Fn key description on 6.5.1 CPXH Series Explanation of Controls . |
| 6 | Set the ultrasonic power level: <ul style="list-style-type: none"> • Press the Fn key until the Power Level icon flashes • Press the Up key to select high power ultrasonics or press the Down key to select low power ultrasonics |
| 7 | Press the Sonics key to activate ultrasonics. |
| 8 | Place the items into a basket, perforated tray, or beakers in a positioning cover. |
| 9 | If using beakers or a solid tray, add cleaning solution to beakers or tray to cover the items. |
| 10 | Slowly lower the tray or beakers into the tank. Do not allow items to contact the tank bottom. |
| 11 | When items are clean, slowly remove them from the tank. |
| 12 | Rinse clean items with clean, warm water and dry, if necessary. |

6.5.6 Cleaning Items (Treating Samples) in Constant Sonics Mode


| NOTICE | |
|---|--|
|  | To stop ultrasonics at any time, press the Sonics key. |

Table 6.20 Treating Samples in Constant Sonics Mode

| Step | Action |
|------|---|
| 1 | Turn Main Power switch on. |
| 2 | Press the On/Standby key to turn on the unit. |
| 3 | If necessary, degas the liquid. See 6.5.3 Degassing for instructions. |
| 4 | Change sonics mode: <ul style="list-style-type: none"> • Press the Fn key until the Constant Sonics icon and the Set Time icon flash • Press the Up key to select Constant Sonics Mode |
| 5 | Set the tank temperature: <ul style="list-style-type: none"> • Press the Fn key until the Set Temperature icon flashes • Press the Up/Down keys to alter the setting to the tank temperature you wish to reach • Press the Heat key once to activate heat. The Heat icon appears <p>NOTICE Units can be switched between °F or °C using the Fn key. See Fn key description on Table 6.15.</p> |
| 6 | Set the ultrasonic power level: <ul style="list-style-type: none"> • Press the Fn key until the Power Level icon flashes • Press the Up key to select high power ultrasonics or press the Down key to select low power ultrasonics |
| 7 | Press the Sonics key to activate ultrasonics. |
| 8 | Place the items into a basket, perforated tray, or beakers in a positioning cover. |
| 9 | If using beakers or a solid tray, add cleaning solution to beakers or tray to cover the items. |
| 10 | Slowly lower the tray or beakers into the tank. Do not allow items to contact the tank bottom. |
| 11 | When items are clean, slowly remove them from the tank. |
| 12 | Rinse clean items with clean, warm water and dry, if necessary. |

6.5.7 Cleaning Items (Treating Samples) in Auto Mode


| NOTICE | |
|---|--|
|  | To stop ultrasonics at any time, press the Sonics key. |

Table 6.21 Treating Samples in Auto Mode

| Step | Action |
|------|--|
| 1 | Turn Main Power switch on. |
| 2 | Press the On/Standby key to turn on the unit. |
| 3 | If necessary, degas the liquid. See 6.5.3 Degassing for instructions. |
| 4 | Set the amount of time you wish the items to be cleaned: <ul style="list-style-type: none"> Use Up/Down keys to increase/decrease cycle time (hold for quick increments/decrements) |
| 5 | Set the tank temperature: <ul style="list-style-type: none"> Press the Fn key until the Set Temperature icon flashes Press the Up/Down keys to alter the setting to the tank temperature you wish to reach Press the Heat key once to activate heat. The Heat icon appears Units can be switched between °F or °C using the Fn key. See Fn key description on 6.5.1 CPXH Series Explanation of Controls. |
| 6 | Set the ultrasonic power level: <ul style="list-style-type: none"> Press the Fn key until the Power Level icon flashes Press the Up key to select high power ultrasonics or press the Down key to select low power ultrasonics |
| 7 | Press the Auto key to begin Auto Cycle. Heater will turn on and sonics will start once set temperature is attained. |
| 8 | Place the items into a basket, perforated tray, or beakers in a positioning cover. |
| 9 | If using beakers or a solid tray, add cleaning solution to beakers or tray to cover the items. |
| 10 | Slowly lower the tray or beakers into the tank. Do not allow items to contact the tank bottom. |
| 11 | When items are clean, slowly remove them from the tank. |
| 12 | Rinse clean items with clean, warm water and dry, if necessary. |


6.5.8 CPXH Temperature Calibration

The CPXH unit temperature measurement is factory calibrated. Use the following instructions to perform periodic calibrations:

Table 6.22 CPXH temperature calibration


| Step | Action |
|------|---|
| 1 | The ultrasonic bath liquid may be at room temperature or may be heated-up to a desired operating temperature (e.g. 40° C). |
| 2 | Press the Fn key 4 times until the Current Temperature icon starts blinking. Press the Up key to select °F. Press the Fn key 2 more times to return to the Ready State. |
| 3 | Press the On/Standby key to turn off the unit. |
| 4 | Simultaneously, press both the On/Standby and the Fn key. Only the bottom left digits and the Current Temperature icon should turn on. |
| 5 | Stir the solution for 15 seconds to ensure thermal uniformity. |
| 6 | Wait 2 minutes after turning the unit on before taking measurements. This allows for the display to be properly updated. |
| 7 | Use the Up/Down keys to change the display temperature to match the actual tank temperature. |
| 8 | Press the On/Standby key to end calibration. |

6.5.9 Draining Your Unit

| WARNING | General Warning |
|---|--|
|  | Do not immerse the unit in water. Unplug the unit from the power source. |

Models 1800 and 2800 do not have a drain. To empty, use the indented side of the rim to pour the used solution into a waste disposal unit, rinse the tank thoroughly and refill with new solution.

Figure 6.5 Draining of Units 1800 and 2800





Models 3800, 5800, and 8800 include a drain and valve kit.

Table 6.23 Draining your ultrasonic bath

| Step | Action |
|------|---|
| 1 | Place the unit to allow easy reach of the drain tube into a waste disposal unit. |
| 2 | Remove the thread protecting cap from the end of the unit's drain pipe. This will expose the white teflon sealing tape on the drain pipe's threads. |

Table 6.23 Draining your ultrasonic bath

| Step | Action |
|------|--|
| 3 | <p>Hand tighten the drain valve onto the drain pipe over the white teflon sealing tape. Finish tightening the valve in place using an adjustable or a 21mm wrench. Tighten the valve no more than one full turn when using the wrench until the handle is on top.</p>  <p>CAUTION</p> <p>Over tightening of the valve can cause damage to the ultrasonic tank. Always use teflon sealing tape or a sealing paste designed for use with stainless steel if retightening or refitting of the drain valve is required.</p> |
| 4 | <p>Hand tighten the hose adaptor into the end of the drain valve. Slide the drain tube over the barbed hose adaptor end.</p>  |
| 5 | <p>Close the drain valve by turning the handle perpendicular to the valve body and the unit is ready to fill with solution. To open the valve and drain the tank, turn the handle so that it is in line with the valve body.</p> |

6.5.10 Solution Temperature Measurement

The following instructions provide an accurate method to obtain consistent thermal measurements using a calibrated temperature measurement instrument. These readings can be used for cleaning process control or to verify the accuracy of the CPXH temperature readings.

Table 6.24 Measuring the Solution Temperature

| Step | Action |
|------|---|
| 1 | Ensure that sonics and heaters are off. |
| 2 | Stir the solutions for 15 seconds to ensure thermal uniformity. |
| 3 | For CPXH models, wait 2 minutes after turning the unit on before taking measurements. This allows for the display to be properly updated. |
| 4 | Suspend a thermocouple in the bath without allowing the probe to touch the tank walls. |

6.5.11 Solution

Table 6.25 Solution Usage


| Item | Definition |
|-------------------|---|
| Solution Activity | The amount of visible activity is not necessarily related to optimum cavitation for cleaning. |
| Degassing 1 | Fresh solutions contain many dissolved gases (usually air), which reduce effective ultrasonic action. Although solutions will naturally degas over time, using Degas Mode speeds up the degassing process. Solutions that have been sitting unused for 24 hours or longer have reabsorbed some gases. |
| Degassing 2 | Degassing mode is also used where gas has to be removed from liquids or samples. |
| Heat | Increases the chemical activity of cleaning solutions. |
| Solvents | Never use solvents. Vapors of flammable solutions will collect under the unit, where ignition is possible from electrical components. |
| Surface Tension | It can be reduced by adding solution to the bath. Reduced surface tension will increase cavitation intensity and enhance cleaning. |
| Renewal | Replace cleaning solutions often to increase ultrasonic cleaning activity. Solutions, as with most chemicals, will become depleted over time. Solutions can become contaminated with suspended soil particles which coat the tank bottom, inhibiting ultrasonic activity. |

6.6 Cleaning Methods

There are two methods of cleaning - direct and indirect. Each has advantages and disadvantages. When in doubt, run test samples using both methods to decide which one produces the best results for you.

6.6.1 Direct Cleaning Method

Figure 6.6 Direct Cleaning Method



How it works:


- Fill the tank with warm water and a cleaning solution
- Place the items to be cleaned in a perforated tray and lower them into the tank. You can also suspend items on a wire and then immerse them in the solution

The advantages of this method are:

- The simplicity of operation
- Cleaning effectiveness

6.6.2 Indirect Cleaning Method

Figure 6.7 Indirect Cleaning Method



How it works:

- Fill the tank with warm water and a cleaning solution. The tank can be filled with any amount of diluted solution as long as it reaches the water level line once the items to be cleaned and accessories are placed into the tank
- Pour your solution medium into one or more beakers or into a solid insert tray
- Place the beakers in a beaker positioning cover or a solid insert tray to fit your unit. Beakers should not touch the tank's bottom


The advantages of this method are:

- Removed soil stays in the beaker or tray so you can easily examine, filter or discard it
- You can use one or more solutions at the same time (two completely different cleaning solutions or one beaker or tray with a cleaning solution and one with a rinse solution)
- Cleaning solution in your tank needs to be changed less often

6.6.3 Non-Cleaning Application

This method can be used for sample preparation, which includes degassing liquids, mixing, homogenization, dissolving solids, cell lysing and dispersion of particles.

Figure 6.8 Non Cleaning Application



How it works:

- Fill the tank with water and wetting agent. The tray or beakers can be filled with any amount of solution as long as the water in the tank outside the tray or beakers reaches the water level line
- Place the beakers in a beaker positioning cover, an insert tray or an insert basket to fit your unit, or place the beakers and flasks onto a Branson support rack. Beakers should not touch the tank's bottom

Chapter 7: Maintenance

| | | |
|------------|--|-----------|
| 7.1 | Optimizing Your Ultrasonic Bath | 58 |
| 7.2 | Troubleshooting | 59 |
| 7.3 | Glass Slide Test | 61 |
| 7.4 | Service Centers | 62 |
| 7.5 | Information for Users on Disposal of Equipment. | 65 |

7.1 Optimizing Your Ultrasonic Bath

Follow these recommendations to optimize your Ultrasonic Bath unit.

Table 7.1 Tanks

| Item | Maintenance |
|--------------------|--|
| Cleaning | Check the tank for contamination whenever you change solution. If necessary, remove contaminants with a nonabrasive cloth and water. |
| Emptying | Always unplug the unit before emptying the tank. Empty the solution into a waste disposal unit. |
| Filling | Always unplug the line cord before filling the tank. Fill the tank to the operating level (with beaker/tray in place), using warm tap water. |
| Low Solution Level | Will cause the unit to fail. When you remove heavy or bulky loads from the tank, the solution level may drop below the operating level. If so, be sure to replace lost solution and degas, if necessary, depending on the amount used. |
| Overload | Do not rest any items on the tank bottom. Weight on the tank bottom dampens sound energy and will cause damage to the transducer. Instead, use a tray and/or beaker positioning cover to support all items. Allow at least 1 inch (2.5 cm) between the tank bottom and the beaker or receptacle for adequate cavitation. |
| Covers | Allow the unit to heat up faster, to a higher temperature, and avoid excessive liquid evaporation. However, obstructing the cover vents will cause the unit to overheat. |

7.2 Troubleshooting

If your unit does not operate satisfactorily, please check the tables below for possible causes before calling your authorized service center.


| WARNING | General Warning |
|---|--|
|  | <p>High voltage inside - dangerous shock hazard.</p> <p>DO NOT attempt to disassemble or repair the unit.</p> |

Table 7.2 Troubleshooting

| Problem | Cause | What to do |
|--|--|--|
| Unit will not start | Unit not plugged in properly. | Plug into functioning electrical outlet. |
| | M/MH - Mechanical timer not ON. | Turn timer clockwise. Press power switch ON. |
| | CPX/CPXH - POWER switch not ON. | Call nearest authorized service center. |
| | CPX/CPXH - On/Standby key malfunctioning. Blown fuse. | Call nearest authorized service center. |
| Unit operates but does not heat solution | Heater malfunctions. | Call nearest authorized service center. |
| | MH - HEAT not ON. CPXH - HEAT not set properly. | Turn heat ON. |
| | CPXH - membrane malfunctioning. | Call nearest authorized service center. |
| Clogged drain | Clogged drain. | Call nearest authorized service center. |
| GFI protected outlet trips | Units may cause GFI circuit trips. | Connect unit to an unprotected outlet. |
| Unit operates, but does not reach set temperature | Malfunctioning heater or sensor components. | Call nearest authorized service center. |
| Err on actual temp and alarm Icon is on. Sonics and Degas operate. Auto and Heat are inactive. | Malfunctioning sensor components. | Call nearest authorized service center. |

Table 7.2 Troubleshooting

| Problem | Cause | What to do |
|--|---|--|
| Unit operates but display does not function. | CPX/CPXH - Control board malfunctioning. | Call nearest authorized service center. |
| Unit stops operating and display is blank with only alarm icon on. | Overheat condition. | For CPXH units, if temperature of 75° C is reached, ultrasonics will stop and won't resume until temperature falls back to 69° C. Check solution level. Refer to 4.3 Temperature for information on over temperature protection. |
| Decreased ultrasonic activity. NOTICE Refer to 7.3 Glass Slide Test for cavitation check. | Solution is not degassed. | Make sure that tank was filled with warm tap water plus cleaning solution and has run 5-10 minutes. |
| | Solution is spent. Solution level is incorrect for load. | Change solution. Adjust solution to within 3/8 inch (1 cm) of the tank's operating level line with load. |
| | Tank bottom is covered with soil particles. | Empty, then clean tank with warm water. Wipe with a nonabrasive cloth. |
| | Using deionized water in the tank. | Deionized water does not cavitate as actively as soapy tap water. |

7.3 Glass Slide Test


Check your ultrasonic bath periodically to test the level of activity of the ultrasonic cavitation. Frequency of testing will depend on your use of the unit, however, Branson suggests running this test monthly.

You will need the following equipment:

- Frosted microscope glass slide (1" x 3"), such as Fisherbrand^{®1} #12-550-343 frosted microscope slides, or equivalent;
- No. 2 lead pencil; and
- General purpose household cleaning solution, such as Dawn^{®2} liquid soap


Test procedure:

1. Prepare a fresh solution with general purpose household cleaning solution (concentration 1%) and warm tap water 49° C – 60° C (120° F – 140° F).
2. Fill the tank to within 3/8 inch (1 cm) of the "operating level" line.
3. Turn the ultrasonics on for at least five to ten minutes to allow for degassing.
4. Prepare the glass slide by first wetting the frosted portion with tap water.



5. With the No. 2 pencil, on the frosted portion make an "X" from corner to corner.
6. Immerse the frosted end of the slide into the solution. Hold the slide vertically and center it in the solution.
7. Make sure that model CPX/CPHX models are in Timed or Constant Sonics Mode, not Degas Mode, then turn ultrasonics On.


The ultrasonics will begin immediately to remove the lead from the slide. All lead should be removed within 10 seconds. If your unit passes this test, its ultrasonic cavitation is acceptable.

| NOTICE | |
|---|---|
|  | <p>To ensure consistency from test to test, be sure to repeat test conditions—use the same solution concentration, liquid level, temperature, type of pencil, length of degassing, etc.</p> |

1. Fisherbrand is a registered trademark of Fisher Scientific Company.
2. Dawn is a registered trademark of Procter & Gamble Company U.S.A.

7.4 Service Centers

With normal use, your Ultrasonic Bath should not require servicing. However, if it fails to operate satisfactorily, first try to diagnose the problem by following the suggestions in the Troubleshooting Guide on [7.2 Troubleshooting](#).

| WARNING | General Warning |
|---|--|
|  | <p>You will void the warranty if you disassemble your unit. High voltage inside the unit is dangerous.</p> |

If you find that your unit needs repair, carefully pack and return it to your local distributor. If under warranty, remember to include proof of purchase. Your unit will be shipped by ground service unless you specify otherwise.

Table 7.3 Authorized Service Centers (North America)

| Name | Address | Tel/Fax Number |
|---------|--|-------------------|
| Branson | c/o Zuniga Logistics 12013 Sara Road Killam Industrial Park Laredo, TX. 78045 | Tel: 877-330-0405 |

Table 7.4 Technical Support (North America)

| Name | Address | Tel/Fax Number |
|---------|---------|--|
| Branson | N/A | Tel: 203-796-0355 Tel: 203-796-0551 |

Table 7.5 Authorized Repair Representatives

| Name | Address | Tel/Fax Number |
|---------------------------------|---|---|
| Alpha Omega Electronics Corp | 2821 National Drive Garland, TX 75041 | Tel: 972-271-5571 Tel: 800-540-4967 Fax: 972-840-3668 |
| Paragon Electronics | 6861 SW 196th Ave, Suite 404 Pembroke Pines, FL 33332 | Tel: 954-434-8191 Fax: 954-434-8385 |

Table 7.5 Authorized Repair Representatives

| Name | Address | Tel/Fax Number |
|---------------------|---|--|
| Crystal Electronics | 1251 Gorham Street Unit 2 Newmarket ON L3Y 8Y6 Canada | Tel: 905-953-9129 Fax: 905-953-7965 |

Table 7.6 Authorized Service Center/Technical Support (Europe)

| Name | Address | Tel/Fax Number |
|------------------------|------------------------------------|--|
| Branson Ultrasonics BV | Vlierberg 26A NL-3755 BS Eemnes | Tel: 31-35-60-98111 Fax: 31-35-60-98120 |


Table 7.7 Authorized Service Center/Technical Support (Asia)

| Name | Address | Tel/Fax Number |
|--|---|---|
| Branson Ultrasonics (Shanghai) Co. Ltd. | 758 East Rong Le Dong Lu Song Jiang Industry Zone Shanghai, 201613 PRC, China | Tel: 86-21-3781-0588 Fax: 86-21-5774-5200 |
| Branson Ultrasonics Asia Pacific Co. Ltd. | Flat A, 5/F Pioneer Building 213 Wai Yip Street Kwun Tong, Kowloon, Hong Kong | Tel: 852-2790-3393 Fax: 852-2790-4998 |
| Branson Ultrasonics Div. of Emerson Electric (India) Pvt. Ltd. | Plot No A 145/6 TTC Industrial Area NIDC Kopar, Navi Mumbai-400705 | Tel: 91-22-64598200/ 220 |
| PT. Global Mega Indonesia | Jl. Jababeka III H Blok C 17 ET Kawasan Industri Jababeka Cikarang Bekasi 17530, Indonesia | Tel: 62-21-8983-6825, Tel. 62-21-8983-6826 Fax: 62-21-8983-6824 |
| Branson Ultrasonics Division of Emerson Japan Ltd. | 4-3-14 Okada, Atsugi-Shi Kanagawa 243-0021, Japan | Tel: 81-46-229-0429 Fax: 81-46-229-0262 |
| Branson Korea Co. Ltd. | DangJeong-dong, 506-7, Gunpo-si, Gyeonggi-do, Republic of Korea | Tel: 82-1577-0631 Fax: 82-31-422-9572 |

Table 7.7 Authorized Service Center/Technical Support (Asia)

| Name | Address | Tel/Fax Number |
|--|--|--|
| Branson Ultrasonics Div. of Emerson Elec (M) Sdn Bhd. | No. 20, Jalan Rajawali 3, Puchong Jaya Industrial Park Batu 8, Jalan Puchong, 47170 Puchong, Selangor, Malaysia | Tel: 603-8076-8608 Fax: 603-8076-8302 |
| Branson Ultrasonics (Philippines Rep Office) | Emerson Building, 104 Laguna Blvd. Laguna Technopark Inc. Sta. Rosa, Laguna Philippines, 4026 | Tel: 63-49-502-8863 Fax: 63-49-502-8860 |
| Branson Ultrasonics Div. of Emerson Electric (South Asia) Pte. Ltd. | Branson Ultrasonics Div. of Emerson Electric (South Asia) Pte. Ltd. | Tel: 65-6891-7600 Fax: 65-6873-7882 |
| Branson Ultrasonics (Taiwan) Division of Emerson Electric Taiwan Co. Ltd. | 5F-3, No. 1, Wu-Chiuan First Road Wu-Ku Ind Zone, Hsin- Chuang City Taipei County, Taiwan, 24892 | Tel: 886-2-2298-0828 Fax: 886-2-2298-9985 |
| Emerson Electric (Thailand) Co. Ltd. | 662/39-40 Rama 3 Rd. Bangpongpan, Yannawa Bangkok, Thailand 10120 | Tel: 662-293-0121-7 Fax: 662-293-0129 |

7.5 Information for Users on Disposal of Equipment



This symbol indicates separate collection of waste electrical and electronic equipment in the EU-countries and EEA (European Economic Area).

Please do not dispose the product with the general household waste. Please use the return and collection system in your country for the disposal of this product.

Index

A

accessories 10
alkaline 22
assistance 2

C

caustic 22
cavitation 8
chemicals 23
cleaning 58
covers 58

D

degas 30
direct cleaning 55
disposal 65
dissolved gases 54
dissolving 8

E

emptying 58

F

FCC regulations 14
filling 58
flammable solutions 32

G

ground leakage current 14

H

harmful chemicals 23

I

IET 5
implode 8
implosions 8
indirect cleaning 55

L

load size 6
low solution level 58

M

mineral acids 32
molded dip 9

MOVs 5

N

non-cleaning application 56

O

overheat 58

overload 58

overvoltages 5

P

personal injury 4

probe 54

property damage 4

R

rinsing items 6

S

safety 2

solution level 6

T

temperature 21

test conditions 61

thermal measurements 54

thermal uniformity 54

transducers 9

U

ultrasonic sound 8

V

vents 4

W

warranty 62

water-based solutions 32