

# Rosemount™ Hx338+

Steam Sterilizable and Autoclavable Pre-Pressurized pH Sensor





# 1 Rosemount Hx338+ IQ/OQ Instruction Sheet

## 1.1 Purpose

The purpose of this document is to describe the procedure for an Installation Qualification (IQ) and Operation Qualification (OQ) for Rosemount™ Hx338+ pH sensors.

Follow these instructions to perform the IQ/OQ at a customer site.

## 1.2 Applicable documents

To carry out the IQ/OQ procedures, the following documents must be available:

- Rosemount™ Hx338+ Sensor [Quick Start Guide](#)
- Reference Manual for the pH transmitter used to conduct the IQ/OQ
- IQ/OQ Form ([IQ/OQ form](#))

## 1.3 Installation Qualification (IQ)

### Procedure

1. Visually inspect the transmitter.
  - a) Check to see if there is any obvious physical damage to the sensor.  
If there is no damage, the sensor passes the Installation Qualification. If there is damage, the sensor fails the Installation Qualification.
  - b) Note any visible defects.
2. Verify documentation.
  - a) Is the final test inspection (Declaration of Quality) for the sensor available?  
Attach a copy to the completed IQ report.
  - b) Is the Rosemount™ 338+ Quick Start Guide available? Attach a copy to the IQ report.
3. Test connection of sensor to transmitter.  
With the sensor connected to a pH transmitter according to the reference manual for each, do the following:
  - a) Remove the protective cap on the sensor and rinse in deionized water.
  - b) Immerse the sensor in pH 4 buffer and swirl.
  - c) Monitor the mV value for at least three minutes until stable.  
mV reading must be positive and over 100 mV to pass.
  - d) Record value and note pass or fail.

If the above steps are completed without any problem, the sensor has passed IQ.

**Postrequisites**

Sign and date on the IQ section of the IQ/OQ Form ([IQ/OQ Form](#)).

## 1.4 Operation Qualification (OQ)

This test verifies the slope and offset of the pH sensor are within specification.

**Procedure**

1. Perform a two-point calibration as directed by the transmitter reference manual. Use buffers of pH 7 and pH 4.
2. Make sure to rinse the sensor carefully when moving the sensor from one buffer to another.
3. Make sure that the mV reading has stabilized in one buffer solution before moving on to the second buffer solution.
4. Once the calibration is completed, verify the sensor slope and offset which should be available from the calibration menu of the transmitter.
  - a) Offset must be between -20 mV and +20 mV to pass.
  - b) Slope must be between 50.3 mV/pH unit and 60.0 mV/pH unit to pass (or 85% to 100% on some transmitters).

If the above steps are completed without any problem, the sensor has passed OQ.

**Postrequisites**

Sign and date on the OQ section of the IQ/OQ Form ([IO/OQ Form](#)).

## 1.5 IQ/OQ form

<b>Installation Qualification</b>		
	<b>Model</b>	<b>Serial Number</b>
<b>Sensor</b>		
<b>Transmitter</b>		

<b>Documentation</b>	<b>Y/N</b>	<b>Comments</b>
<b>Sensor Instruction Manual Available</b>		
<b>Transmitter Instruction Manual Available</b>		

<b>Visual Inspection</b>	<b>Y/N</b>	<b>Comments</b>
<b>Any Damage Noted?</b>		

<b>Test Connection</b>	<b>mV Reading</b>	<b>Pass/Fail (&gt;100 mV = Pass, &lt;100 mV = Fail)</b>
<b>mV Reading in pH4 Buffer</b>		
<b>Comments</b>		

<b>IQ Completion</b>	<b>Y/N</b>	<b>Comments</b>
<b>Passed?</b>		

<b>Signed</b>	
<b>Name</b>	
<b>Title</b>	
<b>Date</b>	

Operation Qualification		
	Model	Serial Number
Sensor		
Transmitter		

Buffers	Part Number
pH 4 Buffer	
pH 7 Buffer	

	Recorded Value	Pass/Fail (Must be between $\pm 20$ mV)
Sensor Offset		

Comments		
----------	--	--

	Recorded Value	Pass/Fail (Must be between 50.3 mV and 59.2 mV per pH unit, or 85% to 100% on some transmitters)
Sensor Slope		

Comments		
----------	--	--

Signed	
Name	
Title	
Date	



#### GLOBAL HEADQUARTERS

Emerson Automation Solutions  
6021 Innovation Blvd  
Shakopee, MN 55379, USA

📞 +1 800 999 9307 or +1 952 906 8888

📠 F +1 952 949 7001

✉️ liquid.csc@emerson.com

#### NORTH AMERICA

Emerson Automation Solutions  
8200 Market Blvd  
Chanhassen, MN 55317

📞 Toll Free +1 800 999 9307

📠 F +1 952 949 7001

✉️ liquid.csc@emerson.com

#### EUROPE

Emerson Automation Solutions  
Neuhofstrasse 19a P.O. Box 1046  
CH-6340 Baar  
Switzerland

📞 T + 41 (0) 41 768 6111

📠 F + 41 (0) 41 768 6300

✉️ liquid.csc@emerson.com

#### MIDDLE EAST AND AFRICA

Emerson Automation Solutions  
Emerson FZE  
Jebel Ali Free Zone  
Dubai, United Arab Emirates, P.O. Box 17033

📞 T +971 4 811 8100

📠 F +971 4 886 5465

✉️ liquid.csc@emerson.com

#### ASIA-PACIFIC


Emerson Automation Solutions  
1 Pandan Crescent  
Singapore 128461  
Singapore

📞 T +65 777 8211

📠 F +65 777 0947

✉️ liquid.csc@emerson.com

 [Linkedin.com/company/Emerson-Automation-Solutions](https://www.linkedin.com/company/Emerson-Automation-Solutions)

 [twitter.com/rosemount\\_news](https://twitter.com/rosemount_news)

 [Facebook.com/Rosemount](https://www.facebook.com/Rosemount)

 [youtube.com/RosemountMeasurement](https://www.youtube.com/RosemountMeasurement)

©2019 Emerson. All rights reserved.

The Emerson logo is a trademark and service mark of Emerson Electric Co. Rosemount is a mark of one of the Emerson family of companies. All other marks are the property of their respective owners.