

LLOYD'S REGISTER TYPE APPROVAL
Issue by: Lloyd's Register Quality Assurance, Inc
Issued to: Emerson Process Management- Valve Automation
For: G Series - G01, G2, G3, G4, G5, G7, G8, G10 and G13

The undernoted documents have been reviewed for compliance with the requirements of the Lloyd's Register Type Approval System Procedure TA14 Version 04 (Sept. 2020) and this Design Appraisal Document forms part of the Certificate.

This certificate is a renewal of certificate 96/60075(E8).

APPROVAL DOCUMENTATION

TEST REPORTS

Type approval application form August 4, 2020
LRQA Inc Houston PRJ11100279303 PQA Assessment 15 Oct 2020
LR Classification Society (China) Co., Ltd. – PQA Assessment QDO2010620 dated 17 Sept 2020



Gopal Nair
Principal Specialist
Lloyd's Register Quality Assurance, Inc.
1330 Enclave Parkway, Suite 200
Houston, Texas 7007

Supplementary Type Approval Terms and Conditions

Type Approval certifies that a representative sample of the product(s) referred to herein has/have been found to meet the applicable design criteria for the use specified herein. It does not mean or imply approval for any other use, nor approval of any product(s) designed or manufactured otherwise than in strict conformity with the said representative sample.

Type Approval is based on the understanding that the manufacturer's recommendations and instructions and any relevant requirements of the Rules and Regulations are complied with.

Type Approval does not eliminate the need for normal inspection and survey procedures required by the Rules and Regulations. Lloyd's Register Quality Assurance Inc. reserves the right to cancel or withdraw this Type Approval Certificate in accordance with the Lloyd's Register Type Approval System Procedure.